

Custom Solution - 8 components 100 µg/mL each of in Methanol

Art. ID LM24-CUS-72322
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bis(2-ethylhexyl) Phthalate-3,4,5,6-d4	[93951-87-2]	100	µg/mL		
	Phthalic acid, benzylbutyl ester (3,4,5,6-D4) (Benzylbutyl phthalate D4)	[93951-88-3]	100	µg/mL		
	Phthalic acid, bis-ethyl ester D4 (Bis-ethyl phthalate-D4, Diethyl phthalate-D4)	[93952-12-6]	100	µg/mL		
	Diisobutyl phthalate-3,4,5,6-d4	[358730-88-8]	100	µg/mL		
	Dimethyl Phthalate-3,4,5,6-d4	[93951-89-4]	100	µg/mL		
	Phthalic acid, bis-butyl ester D4 (Bis-butyl phthalate D4)	[93952-11-5]	100	µg/mL		
	Dicyclohexyl-phthalate D4	[358731-25-6]	100	µg/mL		
	Diethyl phthalate-3,4,5,6-d4	[93952-13-7]	100	µg/mL		