

Custom Solution - 13 components in HNO3 5%

Art. ID LM24-CUS-72326
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	500	mg/L		
	Mercury (Hg)	[7439-97-6]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		