

Custom Solution - 9 components 100 mg/L each of in Acetonitrile

Art. ID LM24-CUS-72340
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis UN1648 II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	100	mg/L		
	Chrysene D12	[1719-03-5]	100	mg/L		
	Benzo(b)fluoranthene D1 2	[93951-98-5]	100	mg/L		
	Naphthalene D8	[1146-65-2]	100	mg/L		
	Perylene D12	[1520-96-3]	100	mg/L		
	Phenanthrene D10	[1517-22-2]	100	mg/L		
	Benzo(a)pyrene D12	[63466-71-7]	100	mg/L		
	Pyrene D10	[1718-52-1]	100	mg/L		
	Indeno(1,2,3-c,d)pyrene D12	[203578-33-0]	100	mg/L		