

**Scopine di(2-thienylglycolate) (2-Thiopheneacetic acid, alpha-hydroxy-alpha-2-thienyl-, (1alpha,2beta,4beta,5alpha,7beta)-9-methyl-3-oxa-9-azatricyclo[3.3.1.0 2,4]non-7-yl ester)**

Art. ID USP-1A04890

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Scopine di(2-thienylglycolate)	[136310-64-0]				