

### **Custom Solution - 4 components in HNO3 2%**

Art. ID LM24-CUS-72363  
Unit 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	500	mg/L		
	Boron (B)	[7440-42-8]	500	mg/L		
	Silicon (Si)	[7440-21-3]	500	mg/L		
	Barium (Ba)	[7440-39-3]	40	mg/L		