

### **Custom Solution - Fe 5 mg/L in HNO3 2%**

Art. ID LM24-CUS-72275-2X100ML  
Unit 2 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### **Description**

Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Iron (Fe)         | [7439-89-6] | 5                   | mg/L |        |        |