

Sorafenib Diarylurea Analog
(4,4'-[carbonylbis(azanediyl)-4,1-phenyleneoxy]]bis(N-methylpyridine-2-carboxamide),
N,N'-Bis{4-[2-(N-methylcarbamoyl)-4-pyridyloxy]phenyl}urea)

Art. ID USP-1A05990

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N,N'-Bis[4-[2-(N-methyl carbamoyl)-4-pyridyloxy]phenyl]urea	[284670-98-0]				