

Low-Energy Charpy V-Notch Specimens (NIST-Verification, 21 °C, 8 mm Striker) - Please ask for availability!

Art. ID	NIST-2561
Unit	set
Deliverydetails	No Dangerous Good /not restricted

Description

The certified value delivered by this Standard Reference Material (SRM) is intended for the indirect verification at 21 °C of Charpy machines equipped with an 8 mm striker, in accordance with the current ASTM Standard E23 and the current ISO Standard 148-2. A unit of SRM 2561 consists of a set of Charpy specimens needed to perform a single indirect verification. SRM 2561 is made from 4340 alloy steel. The bars are finished to length, heat-treated, and machined in SRM specimen lots ranging in size from 1200 to 2000 specimens. Each specimen has a lot number and an identification number (three or four digits).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for current certificate						