

High-Energy Charpy V-Notch Specimens (NIST-Verification, 21 °C, 8 mm Striker) - Please ask for availability!

| | |
|-----------------|-----------------------------------|
| Art. ID | NIST-2563 |
| Unit | set |
| Deliverydetails | No Dangerous Good /not restricted |

Description

The certified value delivered by this Standard Reference Material (SRM) is intended for the indirect verification at 21 °C of Charpy machines equipped with an 8 mm striker, in accordance with the current ASTM Standard E23 and the current ISO Standard 148-2 . A unit of SRM 2563 consists of a set of Charpy specimens needed to perform a single indirect verification. SRM 2563 is made from 4340 alloy steel. The bars are finished to length, heat-treated, and machined in SRM specimen lots ranging in size from 1200 to 2000 specimens. Each specimen has a lot number and an identification number (three or four digits).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------------------------|-------------------|------------|---------------------|------|--------|--------|
| Please ask for current certificate | | | | | | |