

Custom Solution - 16 components in HNO3 5%

Art. ID LM24-CUS-72399
Unit 2 x 50 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cadmium (Cd)	[7440-43-9]	40	mg/L		
	Silver (Ag)	[7440-22-4]	40	mg/L		
	Antimony (Sb)	[7440-36-0]	80	mg/L		
	Arsenic (As)	[7440-38-2]	80	mg/L		
	Nickel (Ni)	[7440-02-0]	80	mg/L		
	Lead (Pb)	[7439-92-1]	80	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	200	mg/L		
	Aluminium (Al)	[7429-90-5]	200	mg/L		
	Iron (Fe)	[7439-89-6]	400	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Boron (B)	[7440-42-8]	1000	mg/L		