

**Custom Solution - 9 components each at 100 mg/L in HNO3 5%, HF 0,5%**

Art. ID LM24-CUS-70721-4X100ML  
Unit 4 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

**Description**

Conservation 24 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		