

Custom Solution - 25 components in HNO3 5%

Art. ID LM24-CUS-71521-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 8 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 20 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 1 | mg/L | | |
| | Boron (B) | [7440-42-8] | 50 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 2 | mg/L | | |
| | Beryllium (Be) | [7440-41-7] | 2 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 0,5 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 2 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 5 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 10 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 20 | mg/L | | |
| | Lithium (Li) | [7439-93-2] | 2,5 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 5 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 2,5 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 2 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 1 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 0,5 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 1 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 2 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 20 | mg/L | | |
| | Tellurium (Te) | [13494-80-9] | 20 | mg/L | | |
| | Titanium (Ti) | [7440-32-6] | 20 | mg/L | | |
| | Thallium (Tl) | [7440-28-0] | 2 | mg/L | | |
| | Uranium (U) | [7440-61-1] | 1 | mg/L | | |
| | Vanadium (V) | [7440-62-2] | 5 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 20 | mg/L | | |