

Custom Solution - 19 components in HNO3 5%

Art. ID LM24-CUS-72426

Unit 1 x 50 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	2000	µg/L		
	Nickel (Ni)	[7440-02-0]	2000	µg/L		
	Arsenic (As)	[7440-38-2]	2000	µg/L		
	Selenium (Se)	[7782-49-2]	200	µg/L		
	Cadmium (Cd)	[7440-43-9]	200	µg/L		
	Tin (Sn)	[7440-31-5]	100000	µg/L		
	Antimony (Sb)	[7440-36-0]	200	µg/L		
	Thallium (Tl)	[7440-28-0]	100	µg/L		
	Lead (Pb)	[7439-92-1]	200	µg/L		
	Uranium (U)	[7440-61-1]	100	µg/L		
	Zinc (Zn)	[7440-66-6]	100000	µg/L		
	Aluminium (Al)	[7429-90-5]	100000	µg/L		
	Phosphorus (P)	[7723-14-0]	2000000	µg/L		
	Manganese (Mn)	[7439-96-5]	100000	µg/L		
	Potassium (K)	[7440-09-7]	4000000	µg/L		
	Copper (Cu)	[7440-50-8]	100000	µg/L		
	Vanadium (V)	[7440-62-2]	100	µg/L		
	Barium (Ba)	[7440-39-3]	100000	µg/L		
	Sodium (Na)	[7440-23-5]	3000000	µg/L		