

Custom Solution - 9 components in HNO3 2%

Art. ID LM24-CUS-72428
Unit 2 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	40	mg/L		
	Iron (Fe)	[7439-89-6]	40	mg/L		
	Aluminium (Al)	[7429-90-5]	40	mg/L		
	Copper (Cu)	[7440-50-8]	80	mg/L		
	Zinc (Zn)	[7440-66-6]	80	mg/L		
	Manganese (Mn)	[7439-96-5]	8	mg/L		
	Chromium (Cr)	[7440-47-3]	4	mg/L		
	Lead (Pb)	[7439-92-1]	4	mg/L		
	Cadmium (Cd)	[7440-43-9]	1	mg/L		