

Custom Solution - 9 components in HNO3 5%

Art. ID LM24-CUS-72429
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle) - 6 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	2000	mg/L		
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	300	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Tin (Sn)	[7440-31-5]	50	mg/L		