

ICP Multi-Element Standard - ICP-MS Tuning Solution - Tune B ICAP - 7 components in 2 %HNO₃

Art. ID LM24-OP-4000084-500ML
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	1	mg/L		
	Bismuth (Bi)	[7440-69-9]	1	mg/L		
	Cerium (Ce)	[7440-45-1]	1	mg/L		
	Cobalt (Co)	[7440-48-4]	1	mg/L		
	Indium (In)	[7440-74-6]	1	mg/L		
	Lithium (Li)	[7439-93-2]	1	mg/L		
	Uranium (U)	[7440-61-1]	1	mg/L		