

**Custom Solution - 23 components in HNO3 1%, HCl 3% + Extended Shelf-Life (Option D):
+ 12 months**

Art. ID LM24-CUS-72447
Unit 2 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months + Extended Shelf-Life (Option D) + 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Arsenic (As)	[7440-38-2]	2	mg/L		
	Boron (B)	[7440-42-8]	20	mg/L		
	Beryllium (Be)	[7440-41-7]	1	mg/L		
	Calcium (Ca)	[7440-70-2]	5000	mg/L		
	Cadmium (Cd)	[7440-43-9]	1	mg/L		
	Cobalt (Co)	[7440-48-4]	1	mg/L		
	Chromium (Cr)	[7440-47-3]	20	mg/L		
	Copper (Cu)	[7440-50-8]	20	mg/L		
	Potassium (K)	[7440-09-7]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Manganese (Mn)	[7439-96-5]	20	mg/L		
	Molybdenum (Mo)	[7439-98-7]	2	mg/L		
	Sodium (Na)	[7440-23-5]	1000	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Phosphorus (P)	[7723-14-0]	500	mg/L		
	Lead (Pb)	[7439-92-1]	20	mg/L		
	Antimony (Sb)	[7440-36-0]	2	mg/L		
	Selenium (Se)	[7782-49-2]	2	mg/L		
	Tin (Sn)	[7440-31-5]	2	mg/L		
	Thallium (Tl)	[7440-28-0]	2	mg/L		
	Vanadium (V)	[7440-62-2]	2	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		