

## Cast Iron, CRM, Block 37 mm x 37 mm x 15 - 20 mm

Art. ID                      CKD-249-Old  
Unit                         block  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,059	%		
	Arsenic (As)	[7440-38-2]	0,018	%		
	Boron (B)	[7440-42-8]	0,018	%		
	Bismuth (Bi)	[7440-69-9]	0,006	%		
	Carbon (C)	[7440-44-0]	3,75	%		
	Cerium (Ce)	[7440-45-1]	0,077	%		
	Cobalt (Co)	[7440-48-4]	0,11	%		
	Chromium (Cr)	[7440-47-3]	0,093	%		
	Copper (Cu)	[7440-50-8]	0,485	%		
	Iron (Fe)	[7439-89-6]	~92,94	%		
	Lanthanum (La)	[7439-91-0]	0,025	%		
	Magnesium (Mg)	[7439-95-4]	0,056	%		
	Manganese (Mn)	[7439-96-5]	0,127	%		
	Molybdenum (Mo)	[7439-98-7]	0,013	%		
	Niobium (Nb)	[7440-03-1]	0,035	%		
	Nickel (Ni)	[7440-02-0]	1,42	%		
	Phosphorus (P)	[7723-14-0]	0,25	%		
	Lead (Pb)	[7439-92-1]	0,023	%		
	Sulfur (S)	[7704-34-9]	0,008	%		
	Antimony (Sb)	[7440-36-0]	0,004	%		
	Selenium (Se)	[7782-49-2]	0,004	%		
	Silicon (Si)	[7440-21-3]	0,34	%		
	Tin (Sn)	[7440-31-5]	0,004	%		
	Titanium (Ti)	[7440-32-6]	0,096	%		
	Vanadium (V)	[7440-62-2]	0,024	%		
	Tungsten (W)	[7440-33-7]	~0,01	%		
	Zinc (Zn)	[7440-66-6]	0,004	%		
	Zirconium (Zr)	[7440-67-7]	0,040	%		