

Custom Solution - 4 components 10 mg/Leach of in HNO3 2%

Art. ID LM24-CUS-72487
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle). (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Titanium (Ti)	[7440-32-6]	10	mg/L		