

**Custom Solution - Ammonium from Ammonium Chloride (NH<sub>4</sub><sup>+</sup> from Ammonium Chloride) 300 mg/L in Water**

Art. ID LM24-CUS-71182-2X100ML  
Unit 2 x 100 mL bottle  
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Ammonium from Ammonium Chloride	Ammonium (NH <sub>4</sub> <sup>+</sup> )	[14798-03-9]	300	mg/L		