

**Trimethoprim Acrylonitrile Analog**  
**(3-(Phenylamino)-2-(3,4,5-trimethoxybenzyl)acrylonitrile)**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | USP-1A05500                       |
| Unit            | 25 mg                             |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

| Text/Information | Analyte/Parameter   | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|---|--------------|---------------------|------|--------|--------|
|                  | (Z)-3-(phenylamino)-2-(<br>3,4,5-trimethoxybenzyl)<br>acrylonitrile | [30078-48-9] |                     |      |        |        |