

Custom Solution - 10 components in 7% HNO₃

Art. ID LM24-CUS-70136-2X500ML
Unit 2 x 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	2	mg/L		
	Iron (Fe)	[7439-89-6]	2	mg/L		
	Manganese (Mn)	[7439-96-5]	2	mg/L		
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Calcium (Ca)	[7440-70-2]	50	mg/L		
	Magnesium (Mg)	[7439-95-4]	50	mg/L		
	Sodium (Na)	[7440-23-5]	50	mg/L		
	Potassium (K)	[7440-09-7]	20	mg/L		
	Strontium (Sr)	[7440-24-6]	1	mg/L		
	Lithium (Li)	[7439-93-2]	1	mg/L		