

**Copper base, Type: Copper, CRM, Rod 3 mm x 100 mm**

Art. ID IMN-CF1-3MM  
Unit rod  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0045	%		
	Arsenic (As)	[7440-38-2]	0,00067	%		
	Bismuth (Bi)	[7440-69-9]	0,00125	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,0042	%		
	Nickel (Ni)	[7440-02-0]	0,00029	%		
	Lead (Pb)	[7439-92-1]	0,0033	%		
	Antimony (Sb)	[7440-36-0]	0,00024	%		
	Tin (Sn)	[7440-31-5]	0,0021	%		
	Zinc (Zn)	[7440-66-6]	0,0057	%		