

Copper base, Type: Copper, CRM, Rod 3 mm x 100 mm

Art. ID IMN-CF6-3MM
Unit rod
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0012	%		
	Arsenic (As)	[7440-38-2]	0,000032	%		
	Bismuth (Bi)	[7440-69-9]	~0,0000012	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,0001	%		
	Nickel (Ni)	[7440-02-0]	~0,00004	%		
	Lead (Pb)	[7439-92-1]	0,00018	%		
	Antimony (Sb)	[7440-36-0]	0,00002	%		
	Tin (Sn)	[7440-31-5]	~0,000006	%		