

**Custom Solution - 2 components 1000 mg/L each of in HNO3 5%, HF 0.1%**

Art. ID LM24-CUS-72519  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

**Description**

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		