

**TAN030/3 - Acid Number (AN) Certified Reference Material, ISO/IEC 17025 & ISO 17034
Certified, nominal value 3.00 mg KOH/g**

Art. ID ARO-TAN030/3
Unit 3 x 50 g
Deliverydetails No Dangerous Good /not restricted

Description

Shelf life: 24 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
nominal value	Potassium Hydroxide	[1310-58-3]	3	mg/g	ASTM D664	