

TBN3 - Base Number (BN) Certified Reference Material, ISO/IEC 17025 & ISO 17034
Certified, nominal value 3.00 mg KOH/g

Art. ID ARO-TBN3
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Description

Shelf life: 24 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
nominal value	Potassium Hydroxide	[1310-58-3]	3	mg/g	ASTM D2896	