

SWS-PMCC10 - Flash Point Secondary Working Standard, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure A, ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-SWS-PMCC10
Unit 3 x 80 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		60	°C	ASTM D93	