

SWS-PMCC20 - Flash Point Secondary Working Standard, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure A, ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-SWS-PMCC20
Unit 3 x 80 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		80	°C	ASTM D93	