

SWS-PMCC20 - Flash Point Secondary Working Standard, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure A, ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-SWS-PMCC20
Unit 3 x 80 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|----------|--------|
| | Flash point | | 80 | °C | ASTM D93 | |