

CRM-PMCC40 - Flash Point Certified Reference Material, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure A, ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-CRM-PMCC40
Unit 3 x 80 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		230	°C	ASTM D93	