

CRM-PMGO01 - Flash Point Certified Reference Material, Diesel, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure A, ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-CRM-PMGO01
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|----------|--------|
| | Flash point | | 60 | °C | ASTM D93 | |