

CRM-PMLU01 - Flash Point Certified Reference Material, Lubricant, Pensky-Martens Closed Cup (PMCC), ASTM D93, Procedure B (Nominal Value (90 °C to 100 °C), ISO/IEC 17025 & ISO 17034 Certified

Art. ID ARO-CRM-PMLU01
Unit 250 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		105	°C	ASTM D93	