

Sulfur Certified Reference Material, ISO/IEC 17025 & ISO 17034 Certified, Sulfur in diesel (synthetic), 30000 mg/kg, 3.0000 % (30000 ppm)

Art. ID ARO-SUSD30000
Unit 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Sulfur (S) | [7704-34-9] | 3 | % | | |