

Sulfur in diesel (synthetic), ASTM D2622, Calibration Set C: 1.0000 % (10000 mg/kg), 2.0000 % (20000 mg/kg), 3.0000 % (30000 mg/kg), 4.0000 % (40000 mg/kg), 5.0000 % (50000 mg/kg)

Art. ID ARO-SUSD-ASTMD2622C
Unit 5 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D2622	