

Sulfur in diesel (synthetic), ASTM D4294, Calibration Set A: Blank, 0.0005 % (5 mg/kg), 0.0010 % (10 mg/kg), 0.0100 % (100 mg/kg), 0.0250 % (250 mg/kg), 0.0500 % (500 mg/kg), 0.0750 % (750 mg/kg), 0.1000 % (1000 mg/kg)

Art. ID ARO-SUSD-ASTMD4294A
Unit 8 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|------------|--------|
| | Sulfur (S) | [7704-34-9] | | % | ASTM D4294 | |