

**Sulfur in diesel (synthetic), ASTM D4294, Calibration Set B: 0.1000 % (1000 mg/kg), 0.2500 % (2500 mg/kg), 0.5000 % (5000 mg/kg), 1.0000 % (10000 mg/kg)**

Art. ID                      ARO-SUSD-ASTMD4294B  
Unit                         4 x 100 mL  
Deliverydetails            No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D4294	