

**Sulfur in diesel (synthetic), ASTM D7039, Calibration Set B: Blank, 0.0025 % (25 mg/kg), 0.0100 % (100 mg/kg), 0.0500 % (500 mg/kg), 0.1000 % (1000 mg/kg), 0.3000 % (3000 mg/kg)**

Art. ID                      ARO-SUSD-ASTMD7039B  
Unit                         6 x 100 mL  
Deliverydetails            No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D7039	