

Sulfur in diesel (synthetic), ASTM D7220, Calibration Set A: Blank, 0.0005 % (5 mg/kg), 0.0010 % (10 mg/kg), 0.0025 % (25 mg/kg), 0.0050 % (50 mg/kg)

Art. ID ARO-SUSD-ASTMD7220A
Unit 5 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D7220	