

**Sulfur in diesel (synthetic), ASTM D7220, Calibration Set B: 0.0050 % (50 mg/kg), 0.0100 % (100 mg/kg), 0.0250 % (250 mg/kg), 0.0500 % (500 mg/kg), 0.0750 % (750 mg/kg), 0.1000 % (1000 mg/kg)**

Art. ID                      ARO-SUSD-ASTMD7220B  
Unit                         6 x 100 mL  
Deliverydetails            No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D7220	