

**Sulfur in diesel (synthetic), IP 336, Calibration Set B: Blank, 0.1000 % (1000 mg/kg), 0.3000 % (3000 mg/kg), 0.5000 % (5000 mg/kg)**

Art. ID                      ARO-SUSD-IP336B  
Unit                         4 x 100 mL  
Deliverydetails            No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	IP 336	