

Sulfur in diesel (synthetic), IP 336, Calibration Set C: 0.5000 % (5000 mg/kg), 1.0000 % (10000 mg/kg), 1.5000 % (15000 mg/kg), 2.0000 % (20000 mg/kg)

Art. ID ARO-SUSD-IP336C
Unit 4 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	IP 336	