

Sulfur in diesel (synthetic), IP 336, Calibration Set E: Blank, 0.0300 % (300 mg/kg), 0.0500 % (500 mg/kg), 0.0700 % (700 mg/kg), 0.1000 % (1000 mg/kg), 0.3000 % (3000 mg/kg), 0.5000 % (5000 mg/kg)

Art. ID ARO-SUSD-IP336E
Unit 7 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	blank	%	IP 336	