

Sulfur in diesel (synthetic), IP 497, Calibration Set A: Blank, 0.0050 % (50 mg/kg), 0.0100 % (100 mg/kg), 0.0250 % (250 mg/kg), 0.0300 % (300 mg/kg), 0.0500 % (500 mg/kg)

Art. ID ARO-SUSD-IP497A
Unit 6 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Sulfur (S) | [7704-34-9] | | % | IP 497 | |