

Sulfur in heavy mineral oil, ASTM D4294, Calibration Set A: Blank, 0.0005 % (5 mg/kg), 0.0010 % (10 mg/kg), 0.0100 % (100 mg/kg), 0.0250 % (250 mg/kg), 0.0500 % (500 mg/kg), 0.0750 % (750 mg/kg), 0.1000 % (1000 mg/kg)

Art. ID ARO-SUHM-ASTMD4294A
Unit 8 x 100 mL
Deliverydetails No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	ASTM D4294	