

**Sulfur in heavy mineral oil, IP 336, Calibration Set F: 0.5000 % (5000 mg/kg), 1.0000 % (10000 mg/kg), 1.5000 % (15000 mg/kg), 2.0000 % (20000 mg/kg), 3.0000 % (30000 mg/kg), 4.0000 % (40000 mg/kg), 5.0000 % (50000 mg/kg)**

Art. ID                      ARO-SUHM-IP336F  
Unit                         7 x 100 mL  
Deliverydetails            No Dangerous Good

Description

Shelf life: 12 Months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]		%	IP 336	