

## **Synthetic Sea Water for ASTM D665 / IP 135**

Art. ID                      ARO-LBR-SSWS-300  
Unit                        300 mL  
Deliverydetails            No Dangerous Good /not restricted

### Description

Shelf life: 12 Months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method    | Source |
|------------------|-------------------|------------|---------------------|------|-----------|--------|
|                  |                   |            |                     |      | ASTM D665 |        |
|                  |                   |            |                     |      | / IP 135  |        |