

Custom Solution - 4 components 10 mg/L each of in HNO3 5%

Art. ID LM24-CUS-72538
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Tellurium (Te)	[13494-80-9]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		