

Custom Solution - 4 components 10 mg/L each of in HNO3 5%

Art. ID LM24-CUS-72538
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Antimony (Sb) | [7440-36-0] | 10 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 10 | mg/L | | |
| | Tellurium (Te) | [13494-80-9] | 10 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 10 | mg/L | | |