

Custom Solution - 7 components 100 mg/L each of in HCl 5%

Art. ID LM24-CUS-72541
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Osmium (Os)	[7440-04-2]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		